


## 200 MHz Variable Gain Photoreceiver



The image shows model OE-300-IN-03-FST with 1.035"-40 threaded flange and coupler ring.


Features

- Adjustable transimpedance gain from 10<sup>2</sup> to 10<sup>8</sup> V/A
- Wide bandwidth up to 200 MHz
- InGaAs-PIN photodiode covering the 800 to 1700 nm wavelength range
- High dynamic input range up to 10 mW optical power
- Very low noise, NEP down to 52 fW/√Hz
- Switchable low pass filters for minimizing wideband noise
- Threaded 1.035"-40 and unthreaded 25 mm dia. free space input available, compatible with many optical standard accessories
- Full manual and remote control capability

Applications

- All-purpose low-noise photoreceiver (O/E converter) for the MHz range
- Time resolved optical pulse and power measurements
- Laser intensity noise measurements (RIN)
- Optical front-end for oscilloscopes, spectrum analyzers, A/D converters and RF lock-in amplifiers

Block Diagram



BS-OE-300-R1

## 200 MHz Variable Gain Photoreceiver

Available Versions

OE-300-IN-03-FST



Internal threaded coupler ring with 30 mm outer diameter (included)

1.035"-40 threaded flange for free space applications compatible with many optical standard accessories.

(Please note: Using the fiber-adapters PRA-FC and PRA-FSMA is not recommended as the small size of the active area can drastically reduce the coupling efficiency.)

OE-300-IN-03-FS



25 mm dia. unthreaded flange for free space applications compatible with many optical standard accessories.

Related OE-300 Models

See separate datasheets for following models on [www.femto.de](http://www.femto.de):

OE-300-SI-10-FST

Si-PIN, 1 x 1 mm, 400 - 1000 nm  
1.035"-40 threaded flange

OE-300-SI-10-FS

Si-PIN, 1 x 1 mm, 400 - 1000 nm  
25 mm dia. unthreaded flange

OE-300-SI-30-FST

Si-PIN,  $\varnothing$  3 mm, 320 - 1000 nm  
1.035"-40 threaded flange

OE-300-SI-30-FS

Si-PIN,  $\varnothing$  3 mm, 320 - 1000 nm  
25 mm dia. unthreaded flange

OE-300-IN-01-FC

InGaAs-PIN,  $\varnothing$  80  $\mu$ m, 900 - 1700 nm  
FC fiber receptacle only

OE-300-S

customized versions available on request

## 200 MHz Variable Gain Photoreceiver

Available Accessories

PRA-PAP



post adapter plate, easy to mount on FEMTO photoreceiver series OE, FWPR, HCA-S and LCA-S

(picture shows model OE-300-SI-10)



PS-15



power supply, input: 100 - 240 VAC, output: ±15 VDC, +400/-250 mA

LUCI-10



compact digital I/O interface for USB remote control, supports opto-isolation of amplifier signal path from PC USB port, 16 digital outputs, 3 opto-isolated digital inputs, bus-powered operation

Specifications

Gain

Test conditions

$V_s = \pm 15 \text{ V}$ ,  $T_A = 25 \text{ }^\circ\text{C}$ , system impedance =  $50 \text{ } \Omega$

Transimpedance gain  
Gain accuracy

$1 \times 10^2 \dots 1 \times 10^8 \text{ V/A}$   
 $\pm 1 \%$

Frequency Response

Lower cut-off frequency  
Upper cut-off frequency

DC/100 Hz, switchable  
up to 200 MHz (see table below),  
switchable to 1 MHz or 10 MHz

Input

Noise equivalent power (NEP)  
Max. CW saturation power

see table below  
see table below

Detector

Detector  
Active area

InGaAs-PIN photodiode  
300  $\mu\text{m}$  dia.

Spectral response  
Sensitivity R  
Dark current

800 - 1700 nm  
0.95 A/W typ. @ 1550 nm  
0.1 nA typ.

## 200 MHz Variable Gain Photoreceiver

Specifications (continued)

Performance Depending on Gain Setting

Gain setting (low noise) (V/A)	10 <sup>2</sup>	10 <sup>3</sup>	10 <sup>4</sup>	10 <sup>5</sup>	10 <sup>6</sup>	10 <sup>7</sup>
Upper cut-off frequency (-3 dB)	200 MHz	80 MHz	14 MHz	3.5 MHz	1.8 MHz	220 kHz
NEP ( $\sqrt{\text{Hz}}$ , @ 1550 nm)	192 pW	23 pW	1.9 pW	410 fW	152 fW	55 fW
Measured at	20 MHz	8 MHz	1.4 MHz	350 kHz	180 kHz	22 kHz
Integrated input noise (RMS)*	4.8 $\mu\text{W}$	370 nW	23 nW	3.4 nW	0.82 nW	64 pW
CW sat. power (@ 1550 nm)	10 mW	1.0 mW	100 $\mu\text{W}$	10 $\mu\text{W}$	1.0 $\mu\text{W}$	100 nW

Gain setting (high speed) (V/A)	10 <sup>3</sup>	10 <sup>4</sup>	10 <sup>5</sup>	10 <sup>6</sup>	10 <sup>7</sup>	10 <sup>8</sup>
Upper cut-off frequency (-3 dB)	175 MHz	80 MHz	14 MHz	3.5 MHz	1.8 MHz	220 kHz
NEP ( $\sqrt{\text{Hz}}$ , @ 1550 nm)	137 pW	6.8 pW	1.4 pW	360 fW	127 fW	52 fW
Measured at	18 MHz	8 MHz	1.4 MHz	350 kHz	175 kHz	22 kHz
Integrated input noise (RMS)*	2.9 $\mu\text{W}$	270 nW	20 nW	3.3 nW	0.82 nW	64 pW
CW sat. power (@ 1550 nm)	1.0 mW	100 $\mu\text{W}$	10 $\mu\text{W}$	1.0 $\mu\text{W}$	100 nW	10 nW

\* The integrated input noise is measured with a shaded input in the full bandwidth ("FBW") setting (referred to 1550 nm). The measurement bandwidth is 3 x the upper cut-off frequency at the specific gain setting; filter slope is a 1<sup>st</sup> order roll-off.

The input referred peak-peak noise can be calculated from the RMS noise as follows:

$$P_{\text{Input noise peak-to-peak}} = P_{\text{Input noise RMS}} \times 6$$

The output noise is given by:

$$U_{\text{Output noise RMS}} = P_{\text{Input noise RMS}} \times \text{gain} \times R$$


$$U_{\text{Output noise peak-to-peak}} = U_{\text{Output noise RMS}} \times 6 = P_{\text{Input noise RMS}} \times \text{gain} \times R \times 6$$

The integrated noise will be reduced considerably by setting the low pass filter to "1 MHz" or "10 MHz" instead of "FBW". This is especially useful for continuous wave (CW) measurements.

Output	Output voltage range	$\pm 1$ V (@ 50 $\Omega$ load), for linear amplification
	Output impedance	50 $\Omega$ (designed for 50 $\Omega$ load)
	Slew rate	1000 V/ $\mu\text{s}$
	Max. output current	$\pm 40$ mA
	Output offset compensation	adjustable by offset potentiometer and external control voltage, output offset compensation range min. $\pm 100$ mV
Ext. Offset Control	Control voltage range	$\pm 10$ V
	Offset control input impedance	15 k $\Omega$
Indicator LED	Function	overload
Digital Control	Control input voltage range	LOW bit: -0.8 ... +1.2 V, HIGH bit: +2.3 ... +12 V
	Control input current	0 mA @ 0 V, 1.5 mA @ +5 V, 4.5 mA @ +12 V
	Overload output	non active: <0.4 V @ 0 ... -1 mA active: typ. 5 ... 5.1 V @ 0 ... 2 mA
Power Supply	Supply voltage	$\pm 15$ V
	Supply current	+110/-90 mA (depends on operating conditions, recommended power supply capability min $\pm 200$ mA)
	Stabilized power supply output	$\pm 12$ V, max. 20 mA, +5 V, max. 150 mA
Case	Weight	320 g (0.74 lb.)
	Material	AlMg4.5Mn, nickel-plated
Input Flange	Material	1.4305 stainless steel, glass bead blasted (1.035"-40 threaded flange) AlMg4.5Mn, nickel-plated (25 mm dia. unthreaded flange)
Coupler Ring	Material	1.4305 stainless steel, glass bead blasted

200 MHz Variable Gain Photoreceiver

Specifications (continued) DC Monitor Output	Monitor output gain	Mode Low noise High speed	Monitor gain Gain setting divided by -1 Gain setting divided by -10
	Monitor output polarity Monitor output voltage range Monitor output bandwidth Monitor output impedance	inverting ±1 V (@ ≥1 MΩ load) DC ... 1 kHz 1 kΩ (designed for ≥1 MΩ load)	
Temperature Range	Storage temperature Operating temperature	-40 ... +80 °C 0 ... +60 °C	
Absolute Maximum Ratings	Max. CW power (averaged) Digital control input voltage Analog control input voltage Power supply voltage	12 mW -5 V/+16 V relative to digital ground DGND (pin 9) ±15 V relative to analog ground AGND (pin 3) ±20 V	
Connectors	Input	OE-300-IN-03-FST OE-300-IN-03-FS	1.035"-40 threaded flange for free space applications and for use with various types of optical standard accessories 25 mm unthreaded round flange for free space applications For optical FC input model see OE-300-IN-01-FC
	Output	BNC jack (female)	
	Power supply	Lemo® series 1S, 3-pin fixed socket (mating plug type: FFA.1S.303.CLAC52) Pin 1: +15 V Pin 2: -15 V Pin 3: GND	
	Control port	Sub-D 25-pin, female, qual. class 2 Pin 1: +12 V (stabilized power supply output) Pin 2: -12 V (stabilized power supply output) Pin 3: AGND (analog ground for pins 1 - 8) Pin 4: +5 V (stabilized power supply output) Pin 5: digital output: overload (referred to pin 3) Pin 6: DC Monitor output Pin 7: NC (= not connected) Pin 8: output offset control voltage input Pin 9: DGND (ground for digital control pins 10 - 16) Pin 10: digital control input: gain, LSB Pin 11: digital control input: gain Pin 12: digital control input: gain, MSB Pin 13: digital control input: AC/DC Pin 14: digital control input: high speed / low noise Pin 15: upper cut-off frequency limit 10 MHz Pin 16: upper cut-off frequency limit 1 MHz Pin 17 - 25: NC (= not connected)	



## 200 MHz Variable Gain Photoreceiver

Scope of Delivery

OE-300-IN-03, threaded coupler ring ("FST" version only), Lemo® 3-pin connector, datasheet, transport package

Remote Control Operation

**General**  
 Remote control input bits are opto-isolated and connected by a logical OR function to the local switch settings. For remote control set the corresponding local switches to "Remote", "DC", "L" (low noise mode) and "FBW", and select the desired setting via a bit code at the corresponding digital inputs.  
 Mixed operation, e.g. local AC/DC setting and remote controlled gain setting, is also possible.


Gain setting	Low noise	High speed	Pin 12 Pin 14=LOW MSB	Pin 11	Pin 10 Pin 14=HIGH LSB
	Gain (V/A)	Gain (V/A)			
	$10^2$	$10^3$	LOW	LOW	LOW
	$10^3$	$10^4$	LOW	LOW	HIGH
	$10^4$	$10^5$	LOW	HIGH	LOW
	$10^5$	$10^6$	LOW	HIGH	HIGH
	$10^6$	$10^7$	HIGH	LOW	LOW
	$10^7$	$10^8$	HIGH	LOW	HIGH

AC/DC setting	Coupling	Pin 13
DC	LOW	
AC	HIGH	

Low pass filter setting	Upper cut-off freq. limit	Pin 15	Pin 16
full bandwidth		LOW	LOW
10 MHz		HIGH	LOW
1 MHz		LOW	HIGH

High speed / low noise setting	Mode	Pin 14
low noise mode		LOW
high speed mode		HIGH

Spectral Responsivity




## 200 MHz Variable Gain Photoreceiver

Typical Performance Characteristic


Frequency response

$$V_{\text{Supply}} = \pm 15 V_{\text{DC}}; R_{\text{Load}} = 50 \Omega$$


Gain setting:  $L10^2, H10^3$




Gain setting:  $L10^3, H10^4$




Gain setting:  $L10^4, H10^5$




Gain setting:  $L10^5, H10^6$



Gain setting:  $L10^6, H10^7$



Gain setting:  $L10^7, H10^8$




## 200 MHz Variable Gain Photoreceiver


Typical Performance  
Characteristic (continued)

Input noise equivalent power (NEP)


Gain setting  $L10^2, H10^3$




Gain setting  $L10^3, H10^4$




Gain setting:  $L10^4, H10^5$




Gain setting:  $L10^5, H10^6$



Gain setting:  $L10^6, H10^7$



Gain setting:  $L10^7, H10^8$






# 200 MHz Variable Gain Photoreceiver

Typical Performance  
Characteristic (continued)


## Signal pulse response

Gain setting L10<sup>2</sup>




Rise: 1.92 ns Fall: 1.93 ns

Gain setting H10<sup>3</sup>




Rise: 2.28 ns Fall: 2.30 ns

Gain setting L10<sup>3</sup>




Rise: 3.45 ns Fall: 3.49 ns

Gain setting H10<sup>4</sup>




Rise: 3.59 ns Fall: 3.61 ns

Gain setting L10<sup>4</sup>



Rise: 25.98 ns Fall: 26.90 ns

Gain setting H10<sup>5</sup>




Rise: 27.24 ns Fall: 27.11 ns

### 200 MHz Variable Gain Photoreceiver


Typical Performance  
Characteristic (continued)

Gain setting L10<sup>5</sup>




Rise: 82.24 ns Fall: 84.12 ns

Gain setting H10<sup>6</sup>




Rise: 81.92 ns Fall: 84.24 ns

Gain setting L10<sup>6</sup>




Rise: 208.6 ns Fall: 212.5 ns

Gain setting H10<sup>7</sup>




Rise: 208.8 ns Fall: 208.9 ns

Gain setting L10<sup>7</sup>



Rise: 1639.2 ns Fall: 1691.2 ns

Gain setting H10<sup>8</sup>




Rise: 1701.6 ns Fall: 1702.4 ns


200 MHz Variable Gain Photoreceiver

Dimensions

Threaded free space input OE-300-IN-03-FST:



Free space input OE-300-IN-03-FS:



FEMTO Messtechnik GmbH  
 Klosterstr. 64  
 10179 Berlin · Germany  
 Phone: +49 30 280 4711-0  
 Fax: +49 30 280 4711-11  
 Email: info@femto.de  
 www.femto.de

Specifications are subject to change without notice. Information provided herein is believed to be accurate and reliable. However, no responsibility is assumed by FEMTO Messtechnik GmbH for its use, nor for any infringement of patents or other rights of third parties which may result from its use. No license is granted by implication or otherwise under any patent or patent rights of FEMTO Messtechnik GmbH. Product names mentioned may also be trademarks used here for identification purposes only.

© by FEMTO Messtechnik GmbH · Printed in Germany