




# Spectral Microscopy

Most flexible  
HSI Camera



## Assessing hyperspectral microscopy in real-time

The FirefLEYE 185 brings hyperspectral imaging to microscopy, and also enables endoscopy. In life sciences the time-saving advantage of the snapshot technology (thanks to no-scanning) really makes itself known, as saving time saves patients. The camera can also monitor real-time processes, such as samples in petri dishes for example.

The relay lens also allows for a lens change without camera re-calibration. And a close-up lens set allows a macroscopic scale view of a spot size of just a few mm.

# FirefLEYE 185

## Advantages

- 3D hyperspectral snapshot imager (x, y,  $\lambda$ )
- VNIR coverage (450-950 nm)
- 125 spectral bands
- 12-bit / 14-bit version
- 25 Hz max frame rate
- Compatible with microscopes



FirefLEYE 185 mounted on a UAS



Hyperspectral Mapping with the FirefLEYE 185

## TECHNICAL SPECIFICATIONS

### FIREFLEYE V185

Technology	Multipoint Spectrometer
Number of Sensors	2
Wavelength Range	450 – 950 nm
Spectral Bands	125
Spectral Sampling	4 nm
Spectral Resolution (FWHM)	8 nm @ 532 nm
Spatial Resolution	50 x 50 pixel
Spatial Resolution 2 <sup>nd</sup> Sensor	1000 x 1000 pixel
Total Spectra / Image	2500
Total Data Points / Cube	0.3 Million
Data Depths	12 bit / 14 bit
Readout	Global shutter
Max Frame Rate	25 Hz
Integration Time	0.1 – 1000 ms
Field of View (FOV)	30°, 20°, 13°, 7°, lens selectable
Power Consumption	7 W
Data Link	2 GigE
Weight	490 g
Size	200 x 67 x 60 mm

## For all scenarios

**UAS Mapping** - The FirefLEYE 185 was the first-ever light-weight hyperspectral snapshot camera used for aerial mapping from a UAS. The user has a choice of lenses, enabling different fields of view for different tasks.

In lab use the FirefLEYE can be equipped with close-up lenses, allowing a macroscopic scale view with a spot size of only a few mm to cm.

Attaching a relay lens to the FirefLEYE provides for full interchangeability to C-mount lenses. Mount the camera on your microscope or endoscope without the need of an additional calibration.

The latest improvements to the FirefLEYE include upgrading the main sensor to a modern CMOS sensor. The 14-bit version has 2x higher signal-to-noise-ratio and 4x higher dynamic range compared to the S185 (needing 4 times longer integration times).